

## Specification For Electrosurgical Generator With Vessel Sealer

### SPECIFICATION

- Advance Cautery Machine (its multipurpose uses)
- **Mode –**
  - Monopolar Mode-**
  - Bipolar Mode –**
    - Vessel Sealer** ( it can seal vessels up to 7mm & Temp. of jaw should not exceed between 60 to 95 degree C during sealing and jaw should be nano coated. )
    - Under Water Bipolar Resection Mode** (basically used by URO and gynae surgeons for TURP or fibroid removal in uterus.
    - Cardio-ablator mode –** (used for cardiac Surgery )
- **Quality Standards Certification -** (US –FDA, EU-CE)
- **Warranty – 3 to 5 year**
- **AMC/CMC Rate –** (AMC- 2% , CMC-5%)

  
01 Apr 24.